



VetTek: The standard because of what comes standard.

Compare **Universal's VetTek system** with either Summit's InnoVet* or InnoVet Select* system and you'll quickly see that it's an unfair comparison. Whether it's our advanced four mA stations, 28" wide table top for larger patients, or simply the spare collimator bulb we include, VetTek's standard configuration gives you more time-saving and performance-enhancing features in every category. If you're in the market for a new x-ray system, the product comparison chart below will show you what many veterinarians already know — **VetTek sets the standard** for excellence.



VetTek vs. InnoVet*			
Feature	Universal's VetTek	InnoVet*	InnoVet Select*
mA Stations	4 Standard, 50/100S, 300/400L	2 Standard, 100S/300L	2 Standard, 100S/300L, 3rd Optional 30S or 400L
Cassette Storage Box	Standard	Not Offered	Option
SurePrep Auto Collimator Light	Standard	Not Offered	Option
SurePrep Auto Tube Prep	Standard	Not Offered	Option
Built-in Spare Collimator Bulb	Standard	Not Offered	Not Offered
Collimator Measuring Tape	Standard	Not Offered	Not Offered
Electric Locks for Tubestand	Standard	Not Offered	Not Offered
Urine Trap	Standard	Not Offered	Not Offered
Cassette Tray Accepts All Sizes (standard or metric)	Standard	Option	Option
Stationary Table Top Dimensions	28" x 56" (1)	24" x 56"	24" x 53"
2-Way Float Top Table	Option - 28" x 56" Top	Not Offered	Option - 33" x 58"
4-Way Float Top Table	Option - 31" x 68" Top	Not Offered	Not Offered
Integrated HF Generator	Option - 500 mA output	Not Offered	Option - only 300 mA output

(1) 28" top allows comfortable positioning of larger patient

*InnoVet and InnoVet Select are trademarks of Summit Industries, Inc. Information current at the time of printing in January 2002. Rev. 1/02 9000-VTcomp



Serving the Veterinarian Profession for Over 70 Years.

www.universalxray.net



11550 West King Street • Franklin Park, IL 60131 • Phone 800.800.6006 • Fax 847.288.7911



See the VetTek at Western Veterinarian Conference Booth#728