

*High frequency
R/F generator*

G100RF



User friendliness and compactness



G100RF is a compact yet powerful and easy-to-use generator designed for effective application in radio/fluoro systems such as tilting and remote controlled tables.

G100RF finds its perfect application in fluoroscopy rooms equipped with every kind of RF tables, such as Mercury, Vision and Viromatic.

The 100KHz High Frequency technology assures the best picture quality under every operating condition along with reduced patient dose and extended X-ray tube lifetime.

Workflow can be maximized by taking advantage of the menu-driven operation and standard Anatomic programs that make the selection of the exposure parameters quick and user friendly.

A graphic digital display shows all the operating information, warning messages and Anatomic Programs selection. For maximum clarity, Rad and fluoro exposure values are shown on large, dedicated LED displays.

The unit can also be configured with an optional AEC (Automatic Exposure Control) to assure perfect exposure in

every exam condition and high speed starter for high workload conditions. One or two tubes versions are available to fit any room configuration

Self-test and full overload protection guarantee a high level of safety for the patient as well as extended life of the entire system.

Reliability meets the highest standards you can expect from this vital component of your examination room. Calibration and service can take profit from a number of technical functions that can also be accessed through a PC-based software.

Villa Sistemi Medicali: long-standing experience at the service of our customers



A DEL Global Technologies Company



VILLA SISTEMI MEDICALI spa
20090 Buccinasco (MI) Italy - Via delle Azalee, 3
Tel: +39-02-48859.1 - Fax: +39-02-4881.844
Internet: www.villasmi.com

COMPANY WITH QUALITY SYSTEM
CERTIFIED BY DNV
ISO 9001/EN 46001